

Pathflow™ Flu A/B Combi (EN)

Instructions for Use

INTENDED USE

The Pathflow™ Flu A/B Combi is a rapid chromatographic immunoassay for the qualitative detection of influenza A and B antigens in nasopharyngeal swab, throat swab or nasal aspirate specimens. It is intended to aid in the rapid differential diagnosis of influenza A and B viral infections.

SUMMARY

Influenza (commonly known as 'flu') is a highly contagious, acute viral infection of the respiratory tract. It is a communicable disease easily transmitted through the coughing and sneezing of aerosolized droplets containing live virus¹. Influenza outbreaks occur each year during the autumn and winter months. Type A viruses are typically more prevalent than type B viruses and are associated with most serious influenza epidemics, while type B infections are usually milder. The gold standard of laboratory diagnosis is 14-day cell culture with one of a variety of cell lines that can support the growth of influenza virus². Cell culture has limited clinical utility, as results are obtained too late in the clinical course for effective patient intervention. Reverse Transcriptase Polymerase Chain Reaction (RT-PCR) is a newer method that is generally more sensitive than culture with improved detection rates over culture of 2-23%³. However, RT-PCR is expensive, complex and must be performed in specialized laboratories. The Pathflow™ Flu A/B Combi qualitatively detects the presence of Influenza A and/or Influenza B antigen in nasopharyngeal swab or throat swab or nasal aspirate specimens, providing results within 15 minutes. The test uses antibodies specific for Influenza A and Influenza B to selectively detect Influenza A and Influenza B antigen in nasopharyngeal swab, throat swab or nasal aspirate specimens.

PRINCIPLE

The Pathflow™ Flu A/B Combi is a qualitative, lateral flow immunoassay for the detection of Influenza A and Influenza B nucleoproteins in nasopharyngeal swab, throat swab or nasal aspirate specimens. In this test, antibodies specific to the Influenza A and Influenza B nucleoproteins is separately coated on the test line regions of the test cassette. During testing, the extracted specimen reacts with the antibodies to Influenza A and/or Influenza B that are coated onto particles. The mixture migrates up the membrane to react with the antibodies to Influenza A and/or Influenza B on the membrane and generate one or two coloured lines in the test regions. The presence of this coloured line in either or both of the test regions indicates a positive result. To serve as a procedural control, a coloured line will always appear in the control region if the test has performed properly.

REAGENTS

The test cassette contains anti-Influenza A and B particles and anti- Influenza A and B coated on the membrane.

MATERIALS

Provided:

Test Cassettes, Extraction Reagent, Extraction Tubes, Sterile Swabs, Package Insert / Instructions for Use, Workstation, Extraction Tube Tips.

Not Provided:

Timer, Aspiration Device.

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use test if pouch is damaged.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.

NFM-019 Version 1

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- The used test should be discarded according to local regulations.

STORAGE AND STABILITY

- Store as packaged at room temperature or refrigerated (2°C-30°C).
- The test is stable through the expiration date printed on the sealed pouch.
- The test must remain in the sealed pouch until use until use.
- **DO NOT FREEZE.**
- Do not use beyond the expiration date.

SPECIMENS

• Nasopharyngeal Swab Sample

Insert a sterilised swab into a nasal cavity securely from a nostril and collect mucocoeperidermis wiping turbinate several times.

• Throat Swab Sample

Insert a sterilised swab into pharynx and collect mucocoeperidermis mainly wiping flare region of post-pharyngeal wall and palatine tonsil several times and be careful not to make saliva attach to the swab.


• Nasal Aspirate

Connect an aspiration catheter to an aspiration trap that is attached to an aspiration device, insert the catheter to nasal cavity from a nostril, start the aspiration device and then collect nasal aspirate sample. Dip a sterilized swab into the collected nasal aspirate sample and make the specimen cling to the swab.


DIRECTIONS FOR USE / INSTRUCTIONS FOR USE

Allow the test, specimen and extraction buffer to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch.
2. Place the Extraction Tube in the workstation. Hold the extraction reagent bottle upside down vertically. Squeeze the bottle and let the solution drop into the extraction tube freely without touching the edge of the tube. Add 10 drops of solution (Approx. 400ul) to the Extraction Tube. See illustration 1.
3. Place the swab specimen in the Extraction Tube. Rotate the swab for approximately 10 seconds while pressing the head against the inside of the tube to release the antigen in the swab. See illustration 2.
4. Remove the swab while squeezing the swab head against the inside of the Extraction Tube as you remove it to expel as much liquid as possible from the swab. Discard the swab in accordance with your biohazard waste disposal protocol. See illustration 3.
5. Fit the dropper tip on top of the extraction tube. Place the test cassette on a clean and level surface. See illustration 4
6. Add three drops of the solution (approx.120ul) to the sample well and then start the timer. Read the result at 15 minutes. Do not interpret the result after 20 minutes.



MICROGEN
BIOPRODUCTS



INTERPRETATION OF RESULTS

	POSITIVE Influenza A: * Two distinct coloured lines appear. One coloured line should be in the control region (C) and another coloured line should be in the Influenza A region (A). A positive result in the Influenza A region indicates that Influenza A antigen was detected in the sample.
	POSITIVE Influenza B: * Two distinct coloured lines appear. One coloured line should be in the control region (C) and another coloured line should be in the Influenza B region (B). A positive result in the Influenza B region indicates that Influenza B antigen was detected in the sample.
	POSITIVE Influenza A and Influenza B: * Three distinct coloured lines appear. One coloured line should be in the control region (C) and two-coloured line should be in the Influenza A region (A) and Influenza B region (B). A positive result in the Influenza A region and Influenza B region indicates that Influenza A antigen and Influenza B antigen were detected in the sample.

***NOTE:** The intensity of the colour in the test line regions (A or B) will vary based on the amount of Flu A or B antigen present in the sample. So, any shade of colour in the test regions (A or B) should be considered positive.

	NEGATIVE: One coloured line appears in the control region (C). No apparent coloured line appears in the test line regions (A or B).
--	--



INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A red line appearing in the control region (C) is an internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative procedural control. If the test is working properly, the background in the result area should be white to light pink and not interfere with the ability to read the test result.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The Pathflow™ Flu A/B Combi is for professional in vitro diagnostic use only. The test should be used for the detection of Influenza A and/or B virus in nasopharyngeal swab, throat swab or nasal aspirate specimens. Neither the quantitative value nor the rate of increase in Influenza A and/or B virus concentration can be determined by this qualitative test.
- The Pathflow™ Flu A/B Combi will only indicate the presence of Influenza A and/or B virus in the specimen from both viable and non-viable Influenza A and B strains.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- A negative result obtained from this kit should be confirmed by culture. A negative result may be obtained if the concentration of the Influenza A and/or B virus present in the nasopharyngeal swab is not adequate or is below the detectable level of the test.
- Excess blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result.
- The accuracy of the test depends on the quality of the swab sample. False negatives may result from improper sample collection or storage.
- The use of over the counter and prescription nasal sprays at high concentrations can interfere with results, leading to either invalid or incorrect test results.
- A positive result for influenza A and/or B does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.

REACTIVITY WITH HUMAN INFLUENZA STRAIN

The Pathflow™ Flu A/B Combi was tested with the following human influenza strains and a discernible line at appropriate test-line regions was observed:

Influenza A Virus	Influenza B Virus
A/NWS/33 10(H1N1)	B/R5
A/Hong Kong/8/68(H3N2)	B/Russia/69
A/Port Chalmers/1/73(H3N2)	B/Lee/40
A/WWS/33(H1N1)	B/Hong Kong/5/72
A/New Jersey/8/76(HswN1)	
A/Mal/302/54(H1N1)	
A/chicken/Yuyao/2/2006 (H5N1)	
A/swine/Hubei/251/2001 (H9N2)	
A/Duck/Hubei/216/1983(H7N8)	
A/Duck/Hubei/137/1982(H10N4)	
A/Anhui/1/2013 (H7N9)	

CLINICAL SENSITIVITY, SPECIFICITY AND ACCURACY

The Pathflow™ Flu A/B Combi has been evaluated with specimens obtained from the patients. RT-PCR is used as the reference method for the Pathflow™ Flu A/B Combi. Specimens were considered positive if RT-PCR indicated a positive result. Specimens were considered negative if RT-PCR indicated a negative result.

Nasopharyngeal Swab Specimen

Pathflow™ Flu A/B Combi / RT/PCR	Type A			Type B			
	RT-PCR		Total	RT-PCR		Total	
	Positive	Negative		Positive	Negative		
Flu A/B Combi	Positive	100	8	102	85	2	87
	Negative	1	180	181	2	200	202
	Total	101	182	283	87	202	289
	Relative Sensitivity	99.0%		97.7%			
	Relative Specificity	98.9%		99.0%			
	Accuracy	98.9%		98.6%			

Throat Swab Specimen

Pathflow™ Flu A/B Combi / RT/PCR	Type A			Type B			
	RT-PCR		Total	RT-PCR		Total	
	Positive	Negative		Positive	Negative		
Flu A/B Combi	Positive	58	1	59	65	1	66
	Negative	3	150	153	4	162	166
	Total	61	151	212	69	163	232
	Relative Sensitivity	95.1%		94.2%			
	Relative Specificity	99.3%		99.4%			
	Accuracy	98.1%		97.8%			

Nasal Aspirate Specimen

Pathflow™ Flu A/B Combi / RT/PCR	Type A			Type B			
	RT-PCR		Total	RT-PCR		Total	
	Positive	Negative		Positive	Negative		
Flu A/B Combi	Positive	46	2	48	94	1	95
	Negative	0	241	241	2	158	160
	Total	46	243	289	96	159	255
	Relative Sensitivity	100%		97.9%			
	Relative Specificity	99.2%		99.4%			
	Accuracy	99.3%		98.8%			

EXPECT VALUE

The Pathflow™ Flu A/B Combi has been compared with a leading commercial RT-PCR test. The correlation between these two systems is over 97%.

PRECISION INTRA-ASSAY & INTER-ASSAY

Within-run and Between-run precision has been determined by using five specimens of Influenza standard control. Three different lots of the Pathflow™ Flu A/B Combi have been tested using negative, Influenza A weak, Influenza B Weak, Influenza A Strong and Influenza B Strong. Ten replicates of each level were tested each day for 3 consecutive days. The specimens were correctly identified >99% of the time.

SPECIFICITY TESTING WITH VARIOUS VIRAL STRAINS

Description	Test Level
Human adenovirus C	5.62 x 10 ⁵ TCID50/ml
Human adenovirus B	1.58 x 10 ⁴ TCID50/ml
Adenovirus type 10	3.16 x 10 ³ TCID50/ml
Adenovirus type 18	1.58 x 10 ⁴ TCID50/ml
Human coronavirus OC43	2.45 x 10 ⁶ LD50/ml
Coxsackievirus A9	2.65 x 10 ⁴ LD50/ml
	1.58 x 10 ⁵ TCID50/ml
Coxsackievirus B5	1.58 x 10 ⁷ TCID50/ml
Human herpesvirus 5	1.58 x 10 ⁴ TCID50/ml
Echovirus 2	3.16 x 10 ⁵ TCID50/ml
Echovirus 3	1 x 10 ⁴ TCID50/ml
Echovirus 6	3.16 x 10 ⁵ TCID50/ml
Herpes simplex virus 1	1.58 x 10 ⁵ TCID50/ml
Human herpesvirus 2	2.81 x 10 ⁵ TCID50/ml
Human Rhinovirus 2	2.81 x 10 ⁴ TCID50/ml
Human Rhinovirus 14	1.58 x 10 ⁶ TCID50/ml
Human Rhinovirus 16	8.89 x 10 ⁵ TCID50/ml
Measles	1.58 x 10 ⁴ TCID50/ml
Mumps	1.58 x 10 ⁴ TCID50/ml
Sendai virus	8.89 x 10 ⁷ TCID50/ml
Parainfluenza virus 2	1.58 x 10 ⁴ TCID50/ml
Parainfluenza virus 3	1.58 x 10 ⁵ TCID50/ml
Respiratory syncytial virus	8.89 x 10 ⁴ TCID50/ml
Human respiratory syncytial virus	1.58 x 10 ⁵ TCID50/ml
Rubella	2.81 x 10 ⁵ TCID50/ml
Varicella-Zoster	1.58 x 10 ⁵ TCID50/ml

TCID50 = Tissue Culture Infectious Dose is the dilution of virus that under the conditions of the assay can be expected to infect 50% of the culture vessels inoculated.

LD50 = Lethal Dose is the dilution of virus that under the conditions of the assay can be expected to kill 50% of the suckling mice inoculated.

CROSS-REACTIVITY

The following organisms were tested at 1.0x10⁸org/ml/ml and all found to be negative when tested with the Pathflow™ Flu A/B Combi:

Arcanobacterium	Pseudomonas aeruginosa
Candida albicans	Staphylococcus aureus subsp. aureus
Corynebacterium	Staphylococcus epidermidis
Enterococcus faecalis	Staphylococcus saprophylicus
Enterococcus faecium	Streptococcus agalactiae
Escherichia coli	Streptococcus bovis
Haemophilus	Streptococcus dysgalatae / subsp. dysgalatae
Moraxella catarrhalis	Streptococcus oralis formerly Streptococcus
Neisseria gonorrhoeae	Streptococcus pneumoniae
Neisseria lactamica	Streptococcus pyogenes
Nisseria subflava	Streptococcus salivarius
Proteus vulgaris	Streptococcus sp group F.type 2

BIBLIOGRAPHY

- Williams, KM, Jackson MA, Hamilton M. (2002) Rapid Diagnostic Testing for URIs in Children; Impact on Physician Decision Making and Cost. *Infect. Med.* 19(3): 109-111.
- Betts, R.F. 1995. Influenza virus, p. 1546-1567. In G.L. Mandell, R.G. Douglas, Jr. and J.E. Bennett (ed.), *Principle and practice of infectious diseases*, 4th ed. Churchill Livingstone, Inc., New York, N.Y.
- WHO recommendations on the use of rapid testing for influenza diagnosis, World Health Organisation, July 2005.

TABLE OF SYMBOLS

SYMBOL	DEFINITION
	Batch Code/Lot Number
	In Vitro Diagnostic Medical Device
	Catalogue Reference
	Upper and Lower Temperature Limit
	Expiry Date
	Manufacturer and Address Details
	Date of Manufacture
	Please consult Instructions for Use
	Sufficient for
	Do not use if package is damaged
	CE Mark
	Caution

This document applies to the following product codes:

M591CE / M591a / M591b

Microgen Bioproducts Limited, Unit 1, Watchmoor Point, Camberley, Surrey, GU15 3AD, United Kingdom.

www: <http://microgenbioproducts.com/>

Email: enquiries.microgen@novacyt.com

Tel: +44 (0) 1276 600081

Fax: +44 (0) 1276 600151

Part of the Novacyt Group

