


**MICROGEN BIOPRODUCTS LTD**



Metapneumovirus Infected R-Mix Too Culture

# **D<sup>3</sup> Metapneumovirus (MPV) Identification Kit**

For the rapid identification of Metapneumovirus directly from  
specimens and in cell culture

**Rapid, accurate results for critical emergent respiratory virus**

**Protecting Food and Health**

---

MICROGEN BIOPRODUCTS LTD

1 Admiralty Way  
Camberley  
Surrey  
GU15 3DT  
United Kingdom

Phone: +44 (0)1276 600081  
Fax: +44 (0)1276 600151  
E-mail: [sales@microgenbioproducts.com](mailto:sales@microgenbioproducts.com)  
[www.microgenbioproducts.com](http://www.microgenbioproducts.com)



## Proven D<sup>3</sup> Testing Methodology for Absolute Confidence in Results

The D<sup>3</sup> DFA MPV identification kit follows the same testing procedure as the industry leading D<sup>3</sup> DFA respiratory virus screening and identification kit.


### Rapid, One Step 15 Minute Stain

Rapid, direct staining of infected cells in patient specimen or culture resulting in fast and accurate detection of Metapneumovirus culture.

### Complete Solution. Everything needed for testing is contained in the D<sup>3</sup> DFA Metapneumovirus Identification Kit

D<sup>3</sup> Metapneumovirus Identification kit Contains:

- 5-ml of MPV FITC-labelled conjugate
- 40x PBS concentrate
- Mounting fluid
- Five control slides with hMPV strains A1, A2, B1 and B2 in separate wells.



Metapneumovirus Infected R-Mix Too Culture

Ordering information:

Kit	
01-030000	D <sup>3</sup> Metapneumovirus Identification Kit