

## Fast, Reliable and Accurate Tests for the Early Detection of Pregnancy

### Rapid Diagnosis:

- 3 minute test result for urine, 5 minutes for serum/plasma.
- Positive results may be achieved within these times.

### User Friendly:

- 36 month expiry.
- Storage at room temperature.
- One step procedure.

### Highly Accurate Results:

- 99.9% sensitivity, specificity and accuracy.
- Urine only test: 25 mIU/ml sensitivity.
- Urine/serum/plasma test: 10 mIU/ml sensitivity.

### Consistent and Reliable:

- Consistent batch to batch performance.
- In-built controls.

### Ordering Information:

Code	Description	Size
60182	hCG rapid test strip (urine) 25 mIU/ml sensitivity	50 test kit
60186	hCG rapid test cassette (urine/serum/plasma) 10 mIU/ml sensitivity	40 test kit



If you need any further information on any of the above tests please contact Microgen Bioproducts Ltd.