

MICROGEN BIOPRODUCTS LIMITED

1 Admiralty Way
Camberley
Surrey GU15 3DT
United Kingdom

Telephone: + 44 (0)1276 600 081
Fax: + 44 (0)1276 600 151



CERTIFICATE OF ANALYSIS

Description: Path-Chek Hygiene Swab PCS-100 **Lot No:** 47737 **Expiry:** 2016-07

Storage: 4 – 25 degrees C **Do not freeze**

Shipping: This product has data supporting stability tolerance during fluctuations in ambient shipping temperatures

TEST PARAMETERS:

SPECIFICATION:

RESULT:

Sterilisation ISOTRN BAND A

For irradiated products, all irradiation batches are challenged with a viable culture and the results of post irradiation sterility check confirm that sterility was achieved.

Swab

Pre-moistened swab with viscose tip and plastic breakaway shaft supplied with sterile peel pouch

Pass

Signed:

A handwritten signature in blue ink, appearing to read 'P. Hughes', is written over the 'Signed:' label.

Date:

17-Aug-15