



## MICROGEN BIOPRODUCTS

### CERTIFICATE OF CONFORMITY

Name of Product: Micropolyspora faeni Kit  
 Product Code: FSK3  
 Components: Micropolyspora faeni antigens, matched positive and negative controls

1. All components in this kit comply with current Microgen Bioproducts approved specifications.
2. Performance complies with current release specifications.

#### RESULTS

1. Fill volume	FSK3a 2.5 - 3.0ml	Complies
	FSK3b 2.5 - 3.0ml	Complies
	FSK3c-d 0.5 - 0.6ml	Complies
2. Azide identity	Present in all reagents	Complies
3. Function	All antigens authenticated by an independant consultant. Positive control - gives precipitin line(s) with all antigens on double diffusion plate. Negative control - gives no precipitin lines with all kit antigens on double diffusion plates.	Complies

#### KEY DATA SUMMARY

Component	Part Number	Lot	Expiry	Lot, expiry correct SOP 340	Component PMP reviewed SOP 337
Negative control	FSKa	30901	2017-03	Pass	Pass
Positive control	FSK3b	As a	As a	Pass	Pass
20mg/ml M.faeni	FSK3c	As a	As a	Pass	Pass
2mg/ml M.faeni	FSK3d	As a	As a	Pass	Pass
Assembly	FSK3	30961	2017-03	Pass	Pass
Overall Expiry Date:		2017-03			

SIGNATURE: 

DATE: 22 / MAR / 2016.