



MICROGEN BIOPRODUCTS

CERTIFICATE OF CONFORMITY

Name of Product:	SALINE CONCENTRATE		
Product Code:	M7		
Batch Number:	<input type="text" value="02961"/>		
Date of Manufacture:	<input type="text" value="2016-10"/>		
Date of Expiry:	<input type="text" value="2021-10"/>		
Components:	Sodium chloride solution.		
Results:			
1. Volume:	Minimum recoverable volume - 150ml.		Complies
2. Appearance:	Colourless liquid		Complies
3. Chloride identity:	Positive		Complies
DATE OF ISSUE:	QUALITY ASSURANCE DEPARTMENT:		
<input type="text" value="17/04/2016"/>	<input type="text" value="RAT"/>		