



MICROGEN BIOPRODUCTS

CERTIFICATE OF CONFORMITY

Name of Product: Micropolyspora faeni Kit
 Product Code: FSK3
 Components: Micropolyspora faeni antigens, matched positive and negative controls

1. All components in this kit comply with current Microgen Bioproducts approved specifications.
2. Performance complies with current release specifications.

RESULTS

1. Fill volume	FSK3a 2.5 - 3.0ml FSK3b 2.5 - 3.0ml FSK3c-d 0.5 - 0.6ml	Complies Complies Complies
2. Azide identity	Present in all reagents	Complies
3. Function	All antigens authenticated by an independant consultant. Positive control - gives precipitin line(s) with all antigens on double diffusion plate. Negative control - gives no precipitin lines with all kit antigens on double diffusion plates.	Complies

KEY DATA SUMMARY

Component	Part Number	Lot	Expiry	Lot, expiry correct SOP 340	Component PMP reviewed SOP 337
Negative control	FSKa	32801	2018-02	PASS	PASS
Positive control	FSK3b	As a	As a	PASS	PASS
20mg/ml M.faeni	FSK3c	As a	As a	PASS	PASS
2mg/ml M.faeni	FSK3d	As a	As a	PASS	PASS
Assembly	FSK3	32871	2018-02	PASS	PASS
Overall Expiry Date:		2018-02			

SIGNATURE: 

DATE: 16/ MAR / 2017