



FINISHED PACK SUMMARY

BATCH NUMBER 9501

FPS M47fCE
Issue Date: 03/20011

TITLE: MICROSCREEN STREP GROUP F LATEX - CE MARK

MICROGEN BIOPRODUCTS

CERTIFICATE OF CONFORMITY

Name of Product: STREP GROUP F CE
Product Code: M47fCE
Components: One of each: M47fCE, Carton WC2268

- 1. This component complies with current Microgen Bioproducts approved specifications.
- 2. Component performance complies with current release specifications.

KEY DATA SUMMARY

Component	Part Number	Lot	Expiry	Lot, expiry correct SOP 340	Component PMP reviewed SOP 337
F Latex	M47fE	9501	2018-09	Pass	Pass
Assembly	M47fCE	9501	2018-09	Pass	Pass

Overall Expiry Date: 2018-09

Component	Appearance	Fill volume	Sodium Azide	Function
F Latex	Milky white suspension	2.75 ± 0.25 ml	Positive	Agglutinates with group F Strep within 30 seconds

COMPLIES

SIGNATURE: [Signature]
DATE: 03 / Jul / 2017