

MICROGEN BIOPRODUCTS LIMITED

1 Admiralty Way
Camberley
Surrey GU15 3DT
United Kingdom



Telephone: + 44 (0)1276 600 081
Fax: + 44 (0)1276 600 151

CERTIFICATE OF ANALYSIS

Description: Path-Chek Listeria Media PC080 **Lot:** 53553 **Expiry:** 2018-07

Storage: 2 — 8 degrees C **Do not freeze**

Shipping: This product has data supporting stability tolerance during fluctuations in ambient shipping temperatures

TEST PARAMETERS:	SPECIFICATION:	RESULT:
Volume :	3.0 +/- 0.2mL	Pass
pH:	7.2 +/- 0.2	Pass
Sensitivity against <i>Listeria monocytogenes</i> :	<10	Pass
Reaction with <i>E coli</i>	negative	Pass
Sterility:	No growth	Pass

Signed: *Sajana*

Quality Assurance

Date: *21 Aug 17*