

Microgen® Strep-ID

An identification system for *Streptococci* and *Enterococci* from clinical, veterinary, food and environmental sources.

- Substrates selected specifically for *Streptococcus*, *Enterococcus* and related species
- Simple and easy-to-use
- Results in 24 hours
- Results interpreted using Microgen ID software

Microgen® Strep-ID is a biochemical test system which utilises a 12 well microwell test strip together with 3 offstrip tests; hippurate hydrolysis (provided), alpha-haemolysis and beta-haemolysis.

As with all Microgen® biochemical ID systems, Strep-ID is supported by the simplest-to-use and most feature-packed Computer Identification System Software (MID60) available, delivering the



most comprehensive and up-to-date databases which registered users can update free-of-charge from the Microgen® Bioproducts website.

Easy to use
Results in 24 hours
CE Marked

Microgen[®] Strep-ID

Species identified:

Streptococcus spp.

- S. acidominimus
- S. agalactiae
- S. anginosus/intermedius
- S. bovis
- S. canis
- S. constellatus
- S. dysgalactiae subsp. dysgalactiae
- S. dysgalactiae subsp. equisimilis
- S. equi subsp. equi
- S. equi subsp. zooepidemicus
- S. equinus
- S. gordonii biovar 1
- S. gordonii biovar 2
- S. gordonii biovar 3
- S. intermedius/anginosus
- S. iniae
- S. mitis biovar 1
- S. mitis biovar 2

- S. mutans
- S. oralis
- S. parasanguinis
- S. pneumoniae
- S. porcinus
- S. pyogenes
- S. salivarius
- S. sanguinis biovar 1
- S. sanguinis biovar 2
- S. sanguinis biovar 3
- S. sanguinis biovar 4
- S. suis
- S. uberis
- S. vestibularis

Enterococcus spp.

- E. avium
- E. casseliflavus
- E. cecorum
- E. dispar
- E. durans
- E. faecalis

- E. faecium
- E. gallinarum
- E. hirae
- E. mundtii
- E. raffinosus

Gemella spp.

- G. haemolysans

Additional reagents available:

- MID61A Nitrate A reagent
- MID61B Nitrate B reagent
- MID61C VP I reagent
- MID61D VP II reagent
- MID61E TDA reagent
- MID61F Indole Kovacs reagent
- MID61G Oxidase strips
- MID61H Mineral oil
- MID61K PYR reagent

Ordering Information:

Code	Description	Size	Storage
MID-62CE	Microgen [®] Strep-ID	20 tests	2 - 8°C

If you need any further information on any of the above tests please contact Microgen Bioproducts Ltd.

M I C R O G E N
B I O P R O D U C T S

Microgen Bioproducts Ltd, Unit 1, Watchmoor Point, Camberley, Surrey GU15 3AD, United Kingdom
T: +44 (0)1276 600 081 F: +44 (0)1276 600 151 E: sales.camberley@lab21.com
www.microgenbioproducts.com

GMM-105-EN revision 2.

Part of the
NOVACYT
G R O U P