

**FINISHED PACK  
SUMMARY**

Page 3 of 3  
FPS M47aCE  
Issue Date: 10/2017

BATCH NUMBER *10401*

TITLE: MICROGEN STREP GROUP A LATEX - CE MARK

***MICROGEN BIOPRODUCTS***

***CERTIFICATE OF CONFORMITY***

Name of Product: STREP GROUP A CE
Product Code: M47aCE
Components: One of each: M47aCE,

1. This component complies with current Microgen Bioproducts approved specifications.
2. Component performance complies with current release specifications.

**KEY DATA SUMMARY**

Component	Part Number	Lot	Expiry	Lot, expiry correct SOP 340	Component PMP reviewed SOP 337
A Latex	M47aCE	<i>10401</i>	<i>2019-09</i>	<i>Pass</i>	<i>Pass</i>
Assembly	M47aCE	<i>10401</i>	<i>2019-09</i>	<i>Pass</i>	<i>Pass</i>

Overall Expiry Date:	<i>2019-09</i>
----------------------	----------------

Component	Appearance	Fill volume	Sodium Azide	Function
A Latex	Milky white suspension	2.75 ± 0.25 ml	Positive	Agglutinates with group A Strep within 30 seconds

**COMPLIES**

SIGNATURE: <i>Sajana</i>
DATE: <i>05 Apr 2018.</i>