

**FINISHED PACK SUMMARY**

BATCH NUMBER 10701

FPS M47dCE  
Issue Date: 11/2017

TITLE: MICROGEN STREP GROUP D LATEX - CE MARK

**MICROGEN BIOPRODUCTS**  
**CERTIFICATE OF CONFORMITY**

Name of Product: STREP GROUP D CE  
Product Code: M47dCE  
Components: One of each: M47dCE

1. This component complies with current Microgen Bioproducts approved specifications.
2. Component performance complies with current release specifications.

**KEY DATA SUMMARY**

Component	Part Number	Lot	Expiry	Lot, expiry correct SOP 340	Component PMP reviewed SOP 337
D Latex	M47dCE	10701	2020-03	PASS	PASS
Assembly	M47dCE	10701	2020-03	PASS	PASS

Overall Expiry Date: 10701

Component	Appearance	Fill volume	Sodium Azide	Function
D Latex	Milky white suspension	2.75 ± 0.25 ml	Positive	Agglutinates with group D Strep within 30 seconds

**COMPLIES**

SIGNATURE: *[Signature]*  
DATE: 29 NOV 18