

**MICROGEN BIOPRODUCTS**

**CERTIFICATE OF CONFORMITY**

**Name of Product:** Microgen Listeria ID 20 test  
**Product Code:** MID-67  
**Components:** 20 x MID 67C; 20 x MID67B; 1 x MID67A

1. All components in this kit comply with current Microgen Bioproducts approved specifications.

2. Performance of strips with QC test panel complies with current release specifications:

L monocytogenes MBCC305	L. seeligeri MBCC310
L. grayii MBCC307	L. welshimeri MBCC308
L. innocua MBCC306	L. ivanovii MBCC325

Give expected pattern of substrate reactions.

3. Product satisfies final Quality Control examination according to FPS MID-67

Component	Inventory Code	Lot Number	Expiry	Lot-expiry correct SOP340	Component PMP reviewed SOP337
Haemolysin Reagent	MID67A	14501A	2020-01	PASS	PASS
Suspending Broth	MID67B	42302	2020-07	PASS	PASS
Test Strip	MID67C	42301	2020-06	PASS	PASS
Assembly	MID67	42392	2020-01	PASS	PASS
<b>Overall Expiry Date:</b>			2020-01		

SIGNATURE: <i>Sajana</i>
DATE: 08 Feb 19