

MICROGEN BIOPRODUCTS

CERTIFICATE OF CONFORMITY

Name of Product: Microgen Listeria ID 20 test
Product Code: MID-67
Components: 20 x MID 67C; 20 x MID67B; 1 x MID67A

- All components in this kit comply with current Microgen Bioproducts approved specifications.
- Performance of strips with QC test panel complies with current release specifications:

L monocytogenes MBCC305	L. seeligeri MBCC310
L grayii MBCC307	L. welshimeri MBCC308
L. innocua MBCC306	L. ivanovii MBCC325

Give expected pattern of substrate reactions.
- Product satisfies final Quality Control examination according to FPS MID-67

Component	Inventory Code	Lot Number	Expiry	Lot-expiry correct SOP340	Component PMP reviewed SOP337
Haemolysin Reagent	MID67A	14901A	2020-06	PASS	PASS
Suspending Broth	MID67B	44601	2020-11	PASS	PASS
Test Strip	MID67C	44601	2021-01	PASS	PASS
Assembly	MID67	44691	2020-06	PASS	PASS
Overall Expiry Date:			2020-06		

SIGNATURE:	<i>Sajana</i>
DATE:	16 Jul '19