

MICROGEN BIOPRODUCTS LIMITED

Unit 1, Watchmoor Point

Camberley

Surrey GU15 3AD

United Kingdom

Telephone: + 44 (0)1276 600 081

Fax: + 44 (0)1276 600 151

MICROGEN
BIOPRODUCTS

CERTIFICATE OF CONFORMITY

Description: Path-Chek Listeria Media **Lot No:** 61103 **Expiry:** 2021-06
Storage: 2 — 8 ° C **Do not freeze**
Shipping: This product has data supporting stability tolerance during fluctuations in ambient shipping temperatures

| TEST PARAMETERS: | SPECIFICATION: | RESULT: |
|---|----------------|---------|
| Volume | 3.0 +/- 0.2mL | Pass |
| pH | 7.2 +/- 0.2 | Pass |
| Sensitivity against <i>Listeria monocytogenes</i> | < 10 | Pass |
| Reaction with <i>E coli</i> | Negative | Pass |
| Sterility | No Growth | Pass |

Signed:

Sanna Hill

Quality Assurance

Date:

28 / Jul / 20