

MICROGEN BIOPRODUCTS LIMITED

Unit 1, Watchmoor Point
Camberley
Surrey GU15 3AD
United Kingdom

MICROGEN
BIOPRODUCTS

Telephone: + 44 (0)1276 600 081
Fax: + 44 (0)1276 600 151

CERTIFICATE OF CONFORMITY

Description: Path-Chek Coliform Media **Lot No:** 62396 **Expiry:** 2022-08
Storage: 2 — 8 ° C **Do not freeze**
Shipping: This product has data supporting stability tolerance during fluctuations in ambient shipping temperatures

| TEST PARAMETERS: | SPECIFICATION: | RESULT: |
|---|----------------|---------|
| Volume | 3.0 +/- 0.2mL | Pass |
| pH | 6.8 +/- 0.2 | Pass |
| Sensitivity against <i>E.Coli</i> | < 10 | Pass |
| Reaction with <i>Citrobacter freundii</i> | Negative | Pass |
| Reaction with <i>Salmonella tranora</i> | Negative | Pass |
| Sterility | No Growth | Pass |

Signed: 

Quality Assurance

Date: 13 April 2021