

MICROGEN

BIOPRODUCTS

CERTIFICATE OF ANALYSIS

DESCRIPTION	LOT NO.	PRODUCT CODE
Path-Chek Coliform Media	63676	PC010
EXPIRY DATE:		2023-01

Test Parameters	Specifications	Result
Volume	3.0 +/- 0.2mL	Pass
pH	6.8 +/- 0.2	Pass
Sensitivity against E.coli	< 10	Pass
Reaction with <i>Citrobacter freundii</i>	Negative	Pass
Reaction with <i>Salmonella tranoroa</i>	Negative	Pass
Sterility	No Growth	Pass

Storage: Store in a refrigerator between 2- 8°C.
Do not freeze.

**We hereby certify that this material/product has been released as
meeting our acceptance criteria**

Signed: *[Signature]*

Date: 21/ Aug/21

Quality Assurance